A woman with long blonde hair, wearing a white strapless dress, stands in a lush, green forest. She is holding a long, winding garland of green leaves and small purple flowers. The scene is filled with soft, glowing light and numerous small, white, star-like sparkles scattered throughout the air and on the foliage. The background is a dense, misty forest with tall trees.

Bring out
your inner
beauty
with
Teoxane's
Beauty
Booster



LABORATORIES

TEOXANE

GENEVA

The logo features three horizontal, overlapping oval shapes of varying lengths, stacked vertically.

mondeal
AESTHETICS

Everybody is talking about hyaluronic acid! But how exactly does it work?

Hyaluronic acid (HA), which is naturally present in our body, performs a fundamental role in keeping the skin in good health. This essential element of the dermis retains water like a sponge so as to maintain its hydration.

SKIN UNDERGOES CHANGES WITH THE PASSAGE OF TIME

The main reasons for the ageing of the skin include excess of oxidants (free radicals)* and progressive loss of intrinsic HA.



With age, we can lose up to two-thirds of our hyaluronic acid pool.¹ At 35 years of age the skin's thickness has decreased on average by about 22%.²

As it gets thinner, the skin becomes more fragile, dehydrates and its elasticity diminishes.

Just like an apple that dries, skin that is deprived of its sources of hydration will lose its shine with the passage of time.

*An unhealthy lifestyle associated with other factors, such as an unbalanced diet and cigarette smoking, affects the skin's appearance and quality.

It's not inevitable! Although our skin undergoes changes, we can take action to stimulate it with Teoxane's Beauty Booster!

AN ALTERNATIVE FOR PREVENTION

The hyaluronic acid injections are designed to prevent premature ageing of our skin and to boost the radiance of our face.³

A COMPENSATION FOR THE SKIN

The hyaluronic acid injections compensate for the natural loss of this substance and stimulate the auto-production of collagen⁴ to help:



The efficacy of injectable hyaluronic acid enables us to postpone the need for anti-ageing surgery.⁵

Approximately 80% of facial skin ageing is attributed to UV-exposure.⁶

TEOXANE'S BEAUTY BOOSTER

Redensity [I]

Holding back the appearance of the first wrinkles

Is your skin starting to change? Would you like to keep it replenished, hydrated and beautiful for longer? Protect and stimulate it with Redensity [I].⁷

A UNIQUE AND PATENTED FORMULA

A high concentration of natural and free hyaluronic acid + a combination of antioxidants, amino acids, minerals and vitamins, all of which are naturally present in our skin.

HYDRATION

HA

ARGININE

ANTIOXIDANT PROTECTION

HA

GLUTATHIONE, N-ACETYL-L-CYSTEIN,
ALPHA-LIPOIC ACID



TISSUE RESTRUCTURATION AND CELL REGENERATION

HA

LEUCINE, ISOLEUCINE, VALINE, GLYCINE, LYSINE,
THREONINE, PROLINE, VITAMIN B₆, ZINC, COPPER

Redensity [I] is a non cross-linked HA gel fluid which allows a smooth application by micro-injections at the chosen areas.

TEOSYAL Redensity® [I] is a product of the TEOSYAL® PureSense range (medical devices, class III, CE mark) for medical professional use only.

Smoothing out a mature skin

Is your skin more creased and thinner than before? Would you like to extend and perfect the filling effect? Hydrate and redensify your skin with Redensity [I].

WHEN AND WHERE SHOULD IT BE USED?⁸


For the face, neck, décolletage and hands. For your personal protocol, ask your practitioner who will recommend the most optimal treatment for your needs.

- Regularly as core treatment to prevent cutaneous ageing
- Alternating with or as a complement to other treatments
- On specific chosen occasions to rehydrate and give a new glow to your skin (for example, before/after prolonged sun exposure...)

VISIBLE ACTION ON CUTANEOUS REDENSIFICATION



BEFORE



Immediately after 2 sessions
of Teosyal® PureSense Redensity I
Quantity injected per session: 2ml

For an optimal result, it is recommended to follow the complete protocol to extend the effects of the treatment and obtain long-lasting results:

- for an optimum treatment: 3 sessions at 3-week intervals for immediate and cumulative results
- for a maintenance treatment of the results, as soon as the need arises: 1 session usually 3 times a year

TEOXANE'S BEAUTY BOOSTER

Redensity [I]

Skin quality improvement ...

IMPROVEMENT:

SKIN GLOW
+60%

ELASTICITY
+56%

TEXTURE
+59%

HYDRATION
+61%

«My impression is that
my face is firmer,
my skin replenished...
I feel better about myself»
Anna, 35

... approved by patients*

91.7%

OF PATIENTS WERE SATISFIED
WITH THEIR TREATMENT

91.7%
ARE READY
TO REPEAT
THE TREATMENT

90.3%
WOULD
RECOMMEND
THE TREATMENT

«I instantly recovered a healthy glow,
everyone told me so»
Olga, 53

«The result is natural and this is
what I was looking for»
Lisa, 47

«Redensity [I] is like a super-cream
that plumps the skin from the inside.
The effect is immediate and long-lasting»
Peggy, 42

* Study report Teosyal® PureSense Redensity [I]. Dr Knoll, Dr Vivier, Dr Carre, Dr Micheels. May 2013. Results from the doctor's evaluation of the global aesthetic improvement for the face indication, 9 weeks after the first injection session. The protocol consisted of 3 injection sessions at 3-week intervals. This study showed a very good treatment tolerance: the main reported side-effects were transient rashes for the face and transient rashes, pain, ecchymosis and itching for the neck and neckline injections.

A woman with long blonde hair, wearing a white strapless dress adorned with a green vine and leaves, stands in a lush, magical forest. The scene is filled with soft, glowing light and numerous small, sparkling particles floating in the air. Purple butterflies are scattered throughout the scene, adding to the ethereal atmosphere. The background consists of tall, dark trees and dense green foliage.

Teoxane's Beauty Booster

TEOSYAL PureSense
Redensity [I]

BRING OUT YOUR INNER BEAUTY

Frequently Asked Questions

WHAT HAPPENS DURING THE SESSION?

- A session consists of intra or hypodermic micro-injections given with a fine needle of very small quantities of Redensity [I] along the lines of wrinkles and evenly spread over the surface of the whole face. In this way the Redensity [I]'s HA and nutrients arrive at the heart of the skin.
- The injection may be manual with syringes and needles, but also electronically assisted for more comfort.
- A session lasts about 20 minutes.

WILL I BE ABLE TO START WORKING AGAIN IMMEDIATELY AFTER A SESSION?

A very slight lack of comfort and a few small spots may appear at the injection points*. These inconveniences generally disappear quickly after injection.



ARE THE RESULTS IMMEDIATELY VISIBLE?

The results are seen immediately at the first session. The skin has a better colour, is less creased, better hydrated and becomes more supple. Gradually, from one session to the next, the complexion looks increasingly radiant and the skin becomes firmer.⁷

DOES IT REQUIRE PAIN RELIEF?

The micro-injections are shallow and performed with very fine needles. For better comfort throughout the injection session Redensity [I] contains lidocaine which is a local anaesthetic widely used in the world. For enhanced comfort it is possible to apply an anaesthetising cream before treatment.

HOW LONG DOES THE EFFECT LAST?

Hyaluronic acid is a resorbable molecule; its effects are temporary and depend on several factors: skin type, age, treated area, injection technique...

For an optimum result, an initial protocol of 3 sessions and a maintenance session 2-3 times a year are recommended.

* There are potential side effects related to the implantation of this device, occurring immediately following injection or after a delay. These include, among others: oedema, erythema, redness, papules. The patient should inform the practitioner of any secondary effect other than those described above or any effect persisting beyond one week.

Located in Geneva, Switzerland, TEOXANE Laboratories are specialised in the design and manufacture of hyaluronic acid-based products for aesthetic medicine.

TEOSYAL PureSense Redensity [I] BEAUTY BOOSTER

Teosyal® PureSense Redensity [I] is indicated for the prevention of wrinkles and rehydration of the face, neck and neckline by injection by a medical professional. Its main components are hyaluronic acid (15mg/ml) and lidocaine hydrochloride (3mg/g). For more information, consult your doctor. This is a regulated health product (class III medical device, CE 0086).

References :

1. Labat-Robert J, et al. Aging of the extracellular matrix. *Médecine & Longévité*. 2012 March; 4(1): 3-32.
2. Farage MA, Miller KW, Elsner P, Maibach HI. Characteristics of the Aging Skin. *Adv Wound Care (New Rochelle)*. 2013 Feb;2(1):5-10.
3. Baspeyras M, et al. Clinical and biometrological efficacy of a hyaluronic acid-based mesotherapy product: a randomised controlled study. *Arch Dermatol Res*. 2013 Oct;305(8):673-82.
4. Anderegg U, et al. More than just a filler - the role of hyaluronan for skin homeostasis. *Exp Dermatol*. 2014 May;23(5):295-303.
5. Data from SNDV (Syndicat National des Dermatologues-Vénérologues)
6. Papakonstantinou E, et al. Hyaluronic acid: A key molecule in skin aging. *Dermatoendocrinol*. 2012 Jul 1;4(3):253-8.
7. Data from AFME (Association Française de Médecine morpho-Esthétique et anti-âge)
8. Teosyal® Puresense Redensity [I] instructions for use.

500310/00 October 2015

 **mondeal**
A E S T H E T I C S

www.lookyounger.com.au


LABORATORIES
TEOXANE
GENEVA